Confidential

No. ZT075-S1059

Jul.20, 2022

NIHON DEMPA KOGYO CO.,LTD.

PRODUCTION DEPT. 3

ULTRASONIC PROBE DEVELOPMENT SECT.

MANAGER Yasunobu Hasegawa
ENGINEER Akihiro Nohara
Kenya Misuo

Report on business requirements

Your Spec No. 11575443 NDK Model No. T01771 NDK Spec No. ZT07S-012

1 Objective Verify business requirements.

	cation detai	İls				
TVP Tag	TPRS Tag	Specification	Test Definition	Sample Size	Verification Method	Criterion
	3.9.1.1-1	Sani-Cloth AF3	Duration: 168 hours at 20 °C (At room temperature)	n=1		
	3.9.1.1-2	Super-Cloth HB	Duration: 168 hours at 20 ºC	n=1		
	3.9.1.1-3	Super Sani-Cloth	Duration: 168 hours at 20 °C (At room temperature)	n=1		
	3.9.1.1-4	Enzol	Duration: 168 hours at 20 ºC	n=1	1	
	3.9.1.1-5	Transeptic	(At room temperature) Duration: 168 hours at 20 °C	n=1		
			(At room temperature) Duration: 500 hours at 20 °C			
	3.9.1.1-6	Metrricide	(At room temperature) Duration: 500 hours at 20 °C	n=1		
	3.9.1.1-7	Dis OPA	(At room temperature)	n=1		
	3.9.1.1-8	Gigasept FF (new)	Duration: 500 hours at 20 °C or 336 hours at 31 °C Submersion Level: Past the strain relief	n=1		1)Transducer shall pass safety testing of
	3.9.1.1-9	Disinfectant : Trophon EPR	Total Cycle: 4,000 cycles	n=1		Hipot and Leakage Current Tests. Hipot (4 000 V [AC rms] / 50 Hz 60
	3.9.1.1-10	RESERT OX HLD	Duration: 168 hours at 30 °C (At room temperature)	n=1	Safety Appearance	seconds) : ≤ 2 mA Leakage Current (264 V [AC rms] / 60 Hz) :
	3.9.1.1-11	Metrricide	Temperature: Room temperature (approx. 20 °C) Total Cycle: 20 cycles Soaking Time: 1 hour per one cycle	n=1	Appearance	Leakage Current (264 V [AC rms] / 60 Hz) : ≤ 0.05 mA 2)Mechanical damage to the plastic housing such as cracking etc.
	3.9.1.1-12	Dis OPA	Temperature: Room temperature (approx. 20 °C) Total Cycle: 20 cycles Soaking Time: 1 hour per one cycle	n=1		
11.0.1	3.9.1.1-13	Gigasept FF (new)	Temperature: Room temperature (approx. 20 °C) Total Cycle: 20 cycles Soaking Time: 1 hour per one cycle	n=1		
	3.9.1.1-14	RESERT OX HLD	Temperature: Room temperature (approx. 20 ºC) Total Cycle: 20 cycles Soaking Time: 1 hour per one cycle	n=1		
	3.9.2.1	Temperature Cycle	Temperature : -10 ~ 60°C Humidity : - Cycle : 20 Cycle Dwell time : 1 hour Ramp up time : - Ramp down Time : - Package : - Operation : -	n=3		1)No mechanical damage to the plastic housing such as cracking etc. 2)Average value of sensitivity :difference before and after the test ±30%(+ 2.27dB/-3.09dB) 3)Hipot (4 000 V [AC rms] / 50 Hz 60 seconds): ≤ 2 mA
	3.9.2.2	WHTOL (Wet High Temperature Operating life Test)	Temperature : 60°C Humidity : 90%RH Cycle : N/A Dwell time : 168 hour Ramp up time : 50 min Ramp down Time : 50min Package : Unpacked Operation : Operation	n=3	Electrical characteristics Safety Appearance	
	3.9.2.3	ISTA 2A	Test Sequence 1) Pre-Humidity 2) Compression 3) 1st Vibration 4) Drop 5) 2nd Vibration Package : Packed Operation : non-Operation	n=3		

2 Verification details

TVP	ation detai	Specification	Test	Sample	Verification Method	Criterion
Tag			Definition The Transducer population	Size		
	3.9.2.4	RDT (Reliability Demonstration Test)	shall have a Reliability of at least 80% @ 80% Confidence Interval over 1,200 hours of operation @ 25 °C ambient	n=1		
	3.9.2.5	4HT (High Heat High Humidity Test)	Temperature : 60°C Humidity : 95%RH Cycle : N/A Dwell time : 250 hour Ramp up time : 50 min Ramp down Time : 50min Package : Unpacked Operation : Non-Operation	n=3		
	3.9.2.6	OTHC (Operating Temperature & Humidity Cycle)	Temperature: +20 ~ +40 °C Humidity: 15 ~ 80% RH Cycle: 10 cycle Dwell time: 1 hour Ramp up time: 50 min Ramp down time: 50 min Package: Unpacked Operation: Operation	n=3	Electrical characteristics Safety	1)No mechanical damage to the plastic housing such as cracking etc. 2)Average value of sensitivity :difference before and after the test ±30%(+2.27dB/-3.09dB)
	3.9.2.7	STHC (Storage Temperature & Humidity Cycle)	Temperature: -40 ~ +60°C and then +10 to +60°C →10 ~ +60°C and then +10 to +60°C Humidity: 10 % RH / 95 % RH Number of Cycles: Minimum of 4 cycles at each set of below Conditions* Dwell time: 1 hour Ramp up time: 50 min Ramp down time: 50 min Packaged: Packed Operation: Non-operating * conditions:	n=3	Appearance	3)Hipot (4 000 V [AC rms] / 50 Hz 60 seconds) : ≤ 2 mA
11.0.1	L		(A) 10°C to 60°C at 95% RH & (B) -40°C to 60°C at 10% RH Up to 10,000 ft : Operating			
	3.9.2.8	Atmospheric Pressure	10,000~36,000 ft: Non- Operating Package: Unpacked Operation: Operation and Non-Operation	n=3	Safety Appearance	No mechanical damage to the plastic housing such as cracking etc. Hipot (4 000 V [AC rms] / 50 Hz 60 seconds) : ≤ 2 mA
	3.9.2.9	Drop Test	Height: 1.22 m Drop: 9 drops(3 drops on of 3 main orientations) Package: Unpacked Operation: Non-Operation	n=3	Safety	Hipot (4 000 V [AC rms] / 50 Hz 60 seconds) : ≤ 2 mA
	3.9.2.10	Vibration	Vibration : Random 5 Hz 0.012 G^2/Hz 50 Hz 0.012 G^2/Hz 100 Hz 0.0006 G^2/Hz Horizontal : 27 minutes, Vertical : 54 minutes	n=3	Safety Appearance	No mechanical damage to the plastic housing such as cracking etc. Average value of sensitivity :difference before and after the test ±30%(+ 2.27dB/-3.09dB) Hipot (4 000 V [AC rms] / 50 Hz 60 seconds) : ≤ 2 mA
	3.9.2.11	Transducer Assembly Reliability and Noise Test (Cable Bending)	Sample level: Transducer Angle: 90° in the initial direction, then 180° in the opposite direction, and then returning 90° to the starting position.(Clockwise and counterclockwise) Centerline of the transducer: 75mm from the point where the transducer strain relief ends on the cable Cycle rate: 20cycles / min. Weight load: 450g Weight location: weight shall be applied the end of the transducer Dwell cycle: 250 k cycles Package: Unpacked Operation: Non-operation Sample size: 3ea	n=3	lmage	Evaluation with ultrasonic diagnostic system: There should be no missing images like black lines in the image in the air.

Siemens Healthineers Ltd.

3 Criterion

Items		Measurement Condition	specification	Criterion
Electrical characteristics	Sensitivity	All elements	Average:-60.0±2.5dB Deviation:(max-min)5.6dB Max.	Average value of sensitivity :difference before and after the test $\pm 30\%(\pm 2.2$ dB/- 3.0 dB)
Safety	Hipot	All channels	AC4KV,cutoff current 2mA,1minite	≦2mA
Salety	Leakage Current	All channels	264V(AC)/60Hz≦50μA	≦ 50μA
Appearance	Appearance	Appearance	Window,Housing, bush,Cable	Mechanical damage to the plastic housing such as cracking etc. No oil leaks.

Resu	

TVP Tag	TPRS Tag	Sample name	Verification Date	Used Equipment	Test Engineer	Result
	3.9.1.1-1	ES1-07	Feb.25, 2022		K.Misuo	Pass
	3.9.1.1-2			Withstanding voltage tester: KIKUSUI		Test omitted because it has the same composition as 3.9.1.1-1
	3.9.1.1-3	ES1-08	Feb.25, 2022	OS8650	K.Misuo	Pass
	3.9.1.1-4	ES1-09	Feb.25, 2022	Tester: sanwa PC7000/PC720M	K.Misuo	Pass
	3.9.1.1-5	ES1-07	Mar.8, 2022]	K.Misuo	Pass
	3.9.1.1-6	ES1-08	Mar.25, 2022		K.Misuo	Pass
	3.9.1.1-7	ES1-09	Mar.25, 2022		K.Misuo	Pass
	3.9.1.1-8	21-002	Jun.23, 2022	According to Siemens Report Technical 9VE4 Transducer disinfection compatible		Pass
	3.9.1.1-9	ES1-15	-	UnderTest		-
	3.9.1.1-10	ES1-07	Mar.16, 2022	Withstanding voltage tester: KIKUSUI TOS8650	K.Misuo	Pass
	3.9.1.1-11	ES1-07	Mar.31, 2022		K.Misuo	Pass
	3.9.1.1-12	ES1-08	Mar.31, 2022	Tester: sanwa PC7000/PC720M	K.Misuo	Pass
	3.9.1.1-13	ES1-14	-	Testing at Siemens		-
11.0.1	3.9.1.1-14	ES1-09	Mar.31, 2022	Withstanding voltage tester: KIKUSUI TOS8650 Tester: sanwa PC7000/PC720M	K.Misuo	Pass
	3.9.2.1	ES1-01 ES1-02 ES1-03	Feb.21, 2022	Temperature Cycle chamber: ESPEC TSE-11-A Chamber: ESPEC SH-641	K.Misuo	Pass
	3.9.2.2	ES1-04 ES1-05 ES1-06	Apr.7, 2022	Pulsar : OLYMPUS 5800PR Digital oscilloscope : Tektoronix DPO3032	K.Misuo	Pass
	3.9.2.3	ES1-01 ES1-02 ES1-03	Apr.2, 2022	Withstanding voltage tester: KIKUSUI TOS8650 Tester: sanwa PC7000/PC720M	K.Misuo	Pass
	3.9.2.4	ES1-10	-		Testing at Siemens	
	3.9.2.5	ES1-01 ES1-02 ES1-03	Apr.22, 2022	Temperature Cycle chamber: ESPEC TSE-11-A Chamber: ESPEC SH-641	K.Misuo	Pass
	3.9.2.6	ES1-04 ES1-05 ES1-06	Apr.15, 2022	Pulsar: OLYMPUS 5800PR Digital oscilloscope: Tektoronix DPO3032	K.Misuo	Pass
	3.9.2.7	ES1-07 ES1-08 ES1-09	May.23, 2022	Withstanding voltage tester: KIKUSUI TOS8650 Tester: sanwa PC7000/PC720M	K.Misuo	Pass
	3.9.2.8	22-002 22-007	Apr.28, 2022	According to Siemens Report SGS-R22-0	0871-EN00	Pass
	3.9.2.9	ES1-01 ES1-02 ES1-03	May.18.2022	Withstanding voltage tester: KIKUSUI TOS8650	K.Misuo	Pass
	3.9.2.10			Accoding to the result of 3.9.2.3		Pass
	3.9.2.11	ES1-11	-	Testing at Siemer	ns	-

No abnormality under the above conditions for the items for which the test was completed.

4 Conclusion

Items that have completed the test have passed the criteria. Items under test will be repoted after the test is completed.

separate sheet1. 3.9.1.1–1 Sani–Cloth AF3
separate sheet2. 3.9.1.1–3 SuperSani–Cloth
separate sheet3. 3.9.1.1–4 Enzol
separate sheet4. 3.9.1.1–5 Transeptic
separate sheet5. 3.9.1.1–6 Metrricide
separate sheet6. 3.9.1.1–7 Dis OPA
separate sheet7. 3.9.1.1–10 RESERT OX HLD
separate sheet8. 3.9.1.1–11 Metrricide
separate sheet9. 3.9.1.1–12 Dis OPA
separate sheet9. 3.9.1.1–12 Dis OPA
separate sheet10. 3.9.1.1–14 RESERT OX HLD
separate sheet11. 3.9.2.1 Temperature Cycle
separate sheet11. 3.9.2.2 WHTOL(Wet High Temperature Operating life Test)
separate sheet13. 3.9.2.3 ISTA 2A
separate sheet14. 3.9.2.5 4HT(High Heat High Humidity Test)
separate sheet15. 3.9.2.6 OTHC(Operating Temperature & Humidity Cycle)
separate sheet16. 3.9.2.7 STHC(Storage Temperature & Humidity Cycle)
separate sheet17. 3.9.2.9 Drop Test

Specification Sani-Cloth AF3

Method Wrap the Sani-Cloth AF3 around the probe and leave it for 168 hours.

After that, check the Hipot, Leakage current and Appearance.

,	<u> </u>	1 1
Used Equipment	Serial No.	Calibration Period
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023



Result

Items	Results and judgment	
Hipot	Pass	
Laskama Cumusint	15.70 <i>μ</i> A	
Leakage Current	Pass	
Appearance	Pass	



Specification Super Sani-Cloth

Method Wrap the Super Sani-Cloth around the probe and leave it for 168 hours.

After that, check the Hipot, Leakage current and Appearance.

Used Equipment	Serial No.	Calibration Period
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023



Result

Items	Results and judgment
Hipot	Pass
Laskara Cumant	15.10 <i>μ</i> A
Leakage Current	Pass
Appearance	Pass



Method

Immerse the probe in Enzol and leave it for 168 hours,

After that, check the Hipot, Leakage current and Appearance.

Used Equipment	Serial No.	Calibration Period
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023



Result

Items	Results and judgment
Hipot	Pass
Laskama Cumusint	16.00 μ A
Leakage Current	Pass
Appearance	Pass



Method Immerse the probe in Transeptic and leave it for 168 hours.

After that, check the Hipot, Leakage current and Appearance.

ľ	Used Equipment	Serial No.	Calibration Period
ľ	TOS8650	12262118	Dec 2022
ŀ	PC7000/PC720M	20025100350	Mar 2023
	PC/000/PC/20M	20025100350	Mar.2023



Result

Items	Results and judgment	
Hipot	Pass	
Leakage Current	15.80 <i>μ</i> A	
	Pass	
Appearance	Pass	



Method

Immerse the probe in Metrricide and leave it for 336 hours.

After that, check the Hipot, Leakage current and Appearance.

Used Equipment	Serial No.	Calibration Period
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023



Result

Items	Results and judgment	
Hipot	Pass	
Leakage Current	16.30 μ A	
	Pass	
Appearance	Pass	



Method

Immerse the probe in Dis OPA and leave it for 500 hours.

After that, check the Hipot, Leakage current and Appearance.

ritter didd, effective inped, Estimage same and rippedianes.		
Used Equipment	Serial No.	Calibration Period
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023



結果

Items	Results and judgment
Hipot	Pass
Leakage Current Appearance	16.90 <i>μ</i> A
	Pass
	Pass



Specification RESERT OX HLD

Method Immerse the probe in RESERT OX HLD and leave it for 168 hours. After that, check the Hipot, Leakage current and Appearance.

	<u>,</u> , ,	
Used Equipment	Serial No.	Calibration Period
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023



Result

Items	Results and judgment	
Hipot	Pass	
Leakage Current	16.30 μ A	
	Pass	
Appearance	Pass	



Method Immerse the probe in Metrricide and perform immersion 1h-5min / 20 cycles. After that, check the Hipot, Leakage current and Appearance.

7 11 001 011010, 0110 011	7 1 11 p c c, = c c 11 c c c c c c c c c c c c c c c	
Used Equipment	Serial No.	Calibration Period
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023



Result

Items	Results and judgment	
Hipot	Pass	
Leakage Current	16.40 μ A	
	Pass	
Appearance	Pass	



Method Immerse the probe in Dis OPA and perform immersion 1h-5min / 20cycles. After that, check the Hipot, Leakage current and Appearance.

7 11 001 011010, 0110 011	7 1 11 p c c, = c c 11 c c c c c c c c c c c c c c c	
Used Equipment	Serial No.	Calibration Period
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023



Result

Items	Results and judgment	
Hipot	Pass	
Leakage Current	14.90 μ A	
	Pass	
Appearance	Pass	



TVP Tag 11.0.1 TPRS Tag 3.9.1.1-14

Specification RESERT OX HLD

Method Immerse the probe in RESERT OX HLD and perform immesion 1h-5min / 20cycles. After that, check the Hipot, Leakage current and Appearance.

Used Equipment	Serial No.	Calibration Period
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023



Result

Items	Results and judgment
Hipot	Pass
Laskara Cumant	15.70 <i>μ</i> A
Leakage Current	Pass
Appearance	Pass



Confidential

Separate sheet 11

TVP Tag 11.0.1 TPRS Tag 3.9.2.1 Specification Temperature Cycle

-10°C/60°C, 1h each, 40cycles Method

After that, check the Electrical characteristics, Hipot, Leakage current and Appearance.

Used Equipment	Serial No.	Calibration Period
TSE-11-A	161000912	Nov.2022
5800PR	40081305	Jun.2022
DPO3032	C010119	Mar.2023
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023

Result

The same	Results and judgment		
Items	ES1-01	ES1-02	ES1-03
Electrical			
characteristics	Pass	Pass	Pass
Hipot	Pass	Pass	Pass
Leakage Current	15.00 <i>μ</i> A	15.10 <i>μ</i> A	15.70 <i>μ</i> A
Leakage Gurrent	Pass	Pass	Pass
Appearance	Pass	Pass	Pass



Separate sheet 12

Specification WHTOL(Wet High Temperature Operating life Test)

Method Temperature: +60°C Humidity: 90%

Cycle:


Dwell time: 168 hours (10080min)

Ramp up time: 50min Ramp down Time: 50min

Package: Unpacked, Operation: Operation

After that, check the Electrical characteristics, Hipot, Leakage current and Appearance.

Used Equipment	Serial No.	Calibration Period
SH-641	92014497	Jun.2022
5800PR	40081305	Jun.2022
DPO3032	C010119	Mar.2023
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023



Result

Items	Results and judgment		
Items	ES1-04	ES1-05	ES1-06
Electrical			
characteristics	Pass	Pass	Pass
Hipot	Pass	Pass	Pass
Leakage Current	14.30 μ A	16.10 <i>μ</i> A	15.50 μ A
Leakage Ourrent	Pass	Pass	Pass
Appearance	Pass	Pass	Pass



Separate sheet 13

Specification ISTA 2A Method 1

1)Pre-Humidity 2)Compression

3)1st Vibration

4)Drop

5)2nd Vibration

Package: Packed, Operation: non-Operation

After that, check the Electrical characteristics, Hipot, Leakage current and Appearance.

Used Equipment	Serial No.	Calibration Period
5800PR	40081305	Jun.2022
DPO3032	C010119	Mar.2023
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar 2023





Result

	D 10 11 1		
Items	Results and judgment		
Items	ES1-01	ES1-02	ES1-03
Electrical			
characteristics	Pass	Pass	Pass
Hipot	Pass	Pass	Pass
Leakage Current	15.30 μ A	15.30 <i>μ</i> A	15.40 μ A
Leanage Ourrent	Pass	Pass	Pass
Appearance	Pass	Pass	Pass



Confidential

Separate sheet 14

Specification 4HT(High Heat High Humidity Test)

Method Temperature : +60°C


Humidity: 95% Cycle: – Dwell time: 250h

Dwell time: 250h
Ramp up time: 50min
Ramp down Time: 50min

Package: Unpacked, Operation: Non-Operation

After that, check the Electrical characteristics, Hipot, Leakage current and Appearance.

Used Equipment	Serial No.	Calibration Period
5800PR	40081305	Jun.2022
DPO3032	C010119	Mar.2023
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023





Result

Items	Results and judgment		
items	ES1-01	ES1-02	ES1-03
Electrical			
characteristics	Pass	Pass	Pass
Hipot	Pass	Pass	Pass
Leakage Current	15.10 μ A	15.40 μ A	15.40 μ A
Leanage Ourrent	Pass	Pass	Pass
Appearance	Pass	Pass	Pass



Confidential

Separate sheet 15

Specification OTHC(Operating Temperature & Humidity Cycle)


Method Temperature: +20~+40°C

Humidity:10~80%
Cycle:10 cycle
Dwell time:1h
Ramp up time:50min
Ramp down time:50min

Package: Unpacked, Operation: Operation

After that, check the Electrical characteristics, Hipot, Leakage current and Appearance.

Used Equipment	Serial No.	Calibration Period
SH-641	92014497	Jun.2022
5800PR	40081305	Jun.2022
DPO3032	C010119	Mar.2023
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023



Result

Items	Results and judgment		
items	ES1-04	ES1-05	ES1-06
Electrical			
characteristics	Pass	Pass	Pass
Hipot	Pass	Pass	Pass
Leakage Current	14.10 μ A	15.00 <i>μ</i> A	15.10 <i>μ</i> A
Leakage Gurrent	Pass	Pass	Pass
Appearance	Pass	Pass	Pass



Separate sheet 16

Specification STHC(Storage Temperature & Humidity Cycle)

Method Temperature: -10~+60°C and then +10~ +60°C

Humidity: 10%/95%

Number of Cycles: Minimum of 4 cycles at each set of below Conditions


(A) +10° C to +60° C at 95% RH & (B) -10° C to +60° C at 10% RH

Dwell time: 1h

Ramp up time: 50min,Ramp down time: 50min Packaged: Packed, Operation: Non-operating

After that, check the Electrical characteristics, Hipot, Leakage current, Appearance.

Used Equipment	Serial No.	Calibration Period
5800PR	40081305	Jun.2022
DPO3032	C010119	Mar.2023
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023





Result

Results and judgment		
ES1-07	ES1-08	ES1-09
Pass	Pass	Pass
Pass	Pass	Pass
16.60 μ A	15.10 μ A	15.10 μ A
Pass	Pass	Pass
Pass	Pass	Pass
	Pass Pass 16.60 μ A Pass	Pass Pass Pass Pass Pass Pass 16.60 \(\mu \) A Pass Pass



Separate sheet 17

Method Height: 1.22 m

Drop: 9 drops(3 drops on of 3 main orientations)

Package: Unpacked Operation : Non-Operation

After that, check the Hipot, Leakage current and Appearance.

 ,	1 , 0	
Used Equipment	Serial No.	Calibration Period
TOS8650	12262118	Dec.2022
PC7000/PC720M	20025100350	Mar.2023

Temperature and humidity at the time of test

	· · · · · · · · · · · · · · · · · · ·
Temperature	Humidity
25.1°C	38.4%

Drop direction1

Drop direction2



Drop direction3



Result

	Result and judgement		
Items	ES1-01	ES1-02	ES1-03
Hipot	Pass	Pass	Pass
Leakage Current	15.30 <i>μ</i> A	15.10 <i>μ</i> A	15.70 <i>μ</i> A
Leakage Gurrent	Pass	Pass	Pass
Appearance	Pass	Pass	Pass

